Gebrauchsanweisung Pflegebett DALI LOW ENTRY 24 Volt Instruction Manual Nursing Bed DALI LOW ENTRY 24 Volt

Modellnummer / Model Number 51.0622.24 S

Typ / Type 181242

((

Stand / Updated **2015.08.28**


Bitte merken:

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Please notice:

In the text of these instructions square brackets and bold figures = **[00]** refer to the parts of the nursing bed.

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- 1 Trapeze (not shown)
- 2 Trapeze pole (not shown)
- 3 End board, foot section
- 4 Release levers for side rails (4)
- 5 Lower leg rest
- 6 Thigh rest
- 7 Side rail beams (4)
- 8 Hand switch
- 9 Backrest
- 10 End board, head section
- 11 Castors with wheel locks (4)
- 12 Reclining surface frame
- 13 Motor drives for backrest and thigh rest
- 14 Transformer unit
- 15 Mattress brackets (4)
- 16 Trapeze pole sockets (2)
- **17** Guide bars (4)

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PREFACE

Dear customer,

Thank you for the confidence you have shown in BURMEIER.

The nursing bed DALI LOW ENTRY 24 Volt is factory-checked for electrical safety and function and has left our factory in faultless condition.

Please take time to read this manual carefully.

This applies especially to the assembly phase but also to daily use.

It is also a practical reference guide. Keep the instructions in a handy place.

These instructions provide information for the operator and user regarding the convenient handling and safe operation.

We wish you and your users much success in care-giving and we are convinced that our product will make a positive contribution.

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GENERAL NOTES

These instructions may describe functions and accessories which are not included with your model.

Overview DALI LOW ENTRY 24 Volt Models:

Model Number: 51.0622.24 S Iron reclining surface, 4 parts Wooden side rails, 2 bars

Before putting the nursing bed into operation for the first time:

- Read these instructions carefully in order to avoid damage or incorrect operation.
- Clean and disinfect the nursing bed before first use.

The bed has been designed with a special height adjustment range of 22-62 cm.

It is therefore important to observe the following instructions when using the bed:

- It is important that the legs of a mobile patient lift do not come into contact with the underside of the bed. Always ensure that there is adequate clearance underneath the bed by adjusting the height of the mattress platform before attempting to transfer a patient.
- Ensure that the caregiver can work comfortably within the bed's height adjustment range or provide appropriate instruction in order to minimise the risk of back injury or strain to the caregiver.

According to the "Medizinprodukte-Betreiberverordnung (MPBetreibV) § 2, before using the nursing bed, the user shall ensure that it is in proper working order and free of defects, and be aware of the instruction manual. This applies also to accessories.

The nursing bed DALI LOW ENTRY 24 Volt complies with all requirements of the guidelines 93/42/EWG for medical products. It is classified as a dynamic medical product of the class I in accordance with "Medizinproduktegesetz (German abbreviation: MPG § 13, Medical Products Act)."

The nursing bed DALI LOW ENTRY 24 Volt is checked by an independent examining institute. As with all technical, electrical devices improper handling can lead to damage and/or injury.

Observe your obligations as operator in accordance with medical devices "Medizinprodukte-Betreiberverordnung (German abbreviation: MPBetreibV, Operators of Medical Products Ordinance)" in order to ensure a dependable and safe operation of this medical device without endangering patients, users and third parties.

This manual contains safety guidelines, which must be observed. All persons who work with the nursing bed DALI LOW ENTRY 24 Volt must familiarize themselves with these instructions and follow the safety guidelines.

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invoivea:

PERSONS

In this manual the following persons are involved:

DEFINITIONS OF INVOLVED

Operator

Operators (e.g.: Medical Supply Houses, Dealers and representatives, health insurance companies) are natural persons or legal entities, who utilize the nursing bed DALI LOW ENTRY 24 Volt or authorize its use. The operator is responsible for instructing the user in the proper use of the unit.

User

Users are persons who as a result of training and experience are authorized to assemble, perform maintenance and/or operate the nursing bed DALI LOW ENTRY 24 Volt. Furthermore, the user can recognize and avoid possible dangers, and judge the clinical condition of the patient.

Patien¹

In this manual a patient is an infirm or handicapped person who reclines on the nursing bed.

Specialist staff

Specialist staff are employees of the operator who are authorized as a result of their education and training to deliver, assemble, dismantle and transport the nursing bed. Furthermore, they are familiar with the regulations for cleaning and disinfecting the unit.

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SAFETY INSTRUCTIONS

The nursing bed DALI LOW ENTRY 24 Volt represents the latest technological development at the time of delivery and is checked by an independent examining institute.

■ Do not use the nursing bed DALI LOW ENTRY 24 Volt if it is damaged.

MEANING OF THE SAFETY SYMBOLS

In these instructions the following safety symbols are used:

Warning about injuries to persons

Dangerous voltage. Life Threatening.

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General danger. Injury or life hazard.

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Warning about property damages

Possibility of damage to drive, material or environment.

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Other symbols

Useful tip. For easier operation or better understanding of the unit.

The safety symbol used does not replace the text of the safety note. Therefore read the safety instructions and follow them exactly!

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SAFETY INSTRUCTIONS FOR OPERATORS

- Remember your responsibilities in accordance with the MPBetreibV in order to ensure the safe operation of this unit without endangering patients, users and third parties.
- Before first use, instruct every user in the safe operation of the nursing bed in accordance with these safety instructions, which must be supplied along with the unit.
- According to MPBetreibV 9 the user shall be notified of the location of this instruction manual.
- Call attention to the dangers of improper use of the unit, especially in regards to the electrical drives and side rails.
- The unit should only be operated by properly trained personnel.
- Make sure that representatives are also familiar with the operation of the nursing bed.
- During long term use it is recommended that the nursing bed be inspected for function and visible damage at regular intervals (Recommendation: yearly) (see chapter 6).
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- When attaching other equipment (e.g. compressors for positioning systems) make sure that they are mounted and can function in a secure and safe manner.
- Do not put multiple outlet electrical extensions under the nursing bed. Leaking liquids can cause a fire hazard.
- Pay special attention to:

Securely locating of all wiring, cables, tubes,

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- ▶ See page 10 : Intended use
- ▶ See page 21 : Safety note for the electric drive system

For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Make sure that your personnel follow the safety instructions.

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SAFETY INSTRUCTIONS FOR USERS

- Let the operator instruct you in the safe use of the nursing bed.
- Before each use make sure that the nursing bed is in proper and faultless condition.
- Take care that no obstacles, such as furniture or slanted ceilings interfere with the adjustment functions.
- Pay attention that when using additional electrical components, such as patient lifts, reading lights or compressors for positioning systems that their electrical cords do not get entangled or damaged by the moving parts of the bed.
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- Make sure that when attaching other appliances (e.g. compressors for positioning systems) that a safe mounting and function is guaranteed.
- Do not put multiple electrical outlets under the nursing bed. Leaking liquids can be a fire
- Pay special attention to securely locating all wiring, cables, tubes, etc.

For further information please contact the manufacturer of the accessories or BURMEIER directly.

- Take the nursing bed out of operation if damage or a malfunction is suspected:
- ► Immediately unplug the electrical plug from the wall outlet;
- ► Clearly mark the nursing bed as "DEFECTIVE";
- ▶ Promptly notify the appropriate operator.

A checklist for assessing the proper condition is located in Chapter 6.

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- that the cable cannot fall or be dragged over the floor.
- Check the power cord regularly for damage (abrasions, exposed wires, kinks, pressure marks, etc.) in particular:
- ▶ after every larger mechanical strain (e.g.: rolling over the power cord with the nursing bed or with an equipment cart; after a strong pulling or bending load caused by the bed rolling away with the power cord still attached to the wall outlet.)
- ▶ after relocating before plugging in the power cable.
- ▶ during prolonged operating by the user.
- Do not use the area under the nursing bed as a "parking spot" for any utensils.
- Lower the reclining surface to its lowest position when leaving the patient unattended. This reduces the risk of injury to the patient getting in and out of the nursing bed.
- Keep the hand switch safe from accidentally falling during non-use (Hang it on the hook). Take care that the the cable is not damaged by the moving parts of the nursing bed
- To protect the patient and especially children from unintentionally adjusting the hand switch, place it beyond their reach (e.g. at the foot end) or lock the hand switch when:
- ▶ the patient is not able to control the bed safely or is unable to get out of dangerous positions without help;
- ▶ the patient could be endangered, because of unintentional adjustments to the drive motor;
- the side rails are up:
 (There is a possible danger that limbs could be squeezed when adjusting the back and thigh rest);
- ▶ unsupervised children are in the room with the nursing bed.
- Adjustment functions can only be performed by a trained person!
- The failure of a lifting motor can result in the bed moving into the Trendelenburg position (head end lowered), which could endanger the patient. If this occurs, the bed must be returned to a horizontal position. Lock the lifting function and inform the operator immediately.

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PRODUCT DESCRIPTION

NORMAL OPERATING PROCEDURE

- The nursing bed DALI LOW ENTRY 24 Volt, hereafter called bed, serves as a comfortable solution for positioning and easing the care of infirm, frail people in senior or nursing homes. Additionally it was developed as a comfortable solution for in-home care of infirm, handicapped or frail people. It will support this care.
- The DALI LOW ENTRY 24 Volt care bed is designed for use in hospitals only in rooms of application group 0.
- This bed can be used for care under the instructions of a physician and can serve in diagnosing, treatment and observation of a patient. It is therefore equipped with a locking mechanism for the electrical adjustment functions.
- This bed has no special connector possibilities for a potential equalization. Consider this before combining with additional electrically operated (medical) devices.

Further references about additionally applicable precautions can be found:

- ▶ In the instructions for the additional electrically operated devices (e.g. compressors for positioning systems, infusion pumps, nutrition
- ▶ In the Standard DIN EN 60601-1-1:2002 (Safety of medical electrical systems)
- ▶ In the Standard VDE 0107:1994 (High current systems in hospitals)
- The maximum long-term load for this bed is 185 kg (patient and accessories). Under this load, the bed can be operated without restric-
- The permitted weight of the patient depends on the total weight of the accessories attached at the same time (e.g. respirators, infusions,...)

Example:

Weight of accessories (with mattress):	Max. permitted weight of patient				
10 kg	165 kg				
40 kg	135 kg				

- ▶ This bed is not appropriate for patients less than 150 cm tall. Observe, especially with patients in poor health condition, the safety tips in chapter 4.
- This bed may only be operated by trained personnel.

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- This bed is useful for multiple re-use. Bear in mind the necessary requirements:
- ▶ See chapter 5;
- ▶ See chapter 6.
- While the patient also is lying in bed, it is allowed to drive the bed within the building.

This bed may only be operated under the operating instructions laid out in this instruction manual.

Any other use is deemed to be inappropriate.

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STRUCTURAL DESIGN

The bed is shipped disassembled on a storage aid. It can be transported without problems even in apartment homes. It consists of two end boards (head- and foot section); a reclining surface frame which is divided in the middle; four side rail beams and a trapeze pole with trapeze handle. The bed stands on four steerable, individually-locking casters.

Reclining surface frame

The reclining surface frame has four sections: a movable backrest, a firm middle section, a moveable thigh rest and lower leg rest. The back- and thigh rest can be adjusted with the aid of electric motors. The horizontal height of the reclining surface can be adjusted and inclined to a low foot end position. All adjustment functions are controlled by a hand switch.

Side Rails

The bed has side rail beams on both sides which can be raised as a barrier or lowered when not needed. This safeguards the patient from accidentally falling out of the bed.

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Electric adjustment system

The electrical adjustment system of this bed is error protected, flame-retardant (V0) and consists of:

■ an external transformer unit.

The transformer unit consists of a power cable, the transformer, and a low voltage connection cable.

The transformer creates a 24 volt low voltage, which is safe for patient and user. The transformer supplies all drive motors with the 24 volt safety low voltage through cables and a wiring harness. The connector on the base frame is moisture proof.

■ a hand switch with a sturdy hook.

The user can lock the adjustment functions of the hand switch if the poor clinical condition of the patient warrants is.

- two drive motors with integrated control for the horizontal height .
- one drive motor with integrated control for the thigh rest .
- one drive motor with integrated control for the backrest.

DESCRIPTION OF MATERIALS

The bed is mostly constructed out of steel sections, their surfaces have been polyester-powder-coated, or with a metallic coating of zinc and chrome. The end boards, the side rail beams and the spring wood slats consist of wood and derived timber products whose surface has been sealed.

All surfaces are harmless to skin contact.

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ASSEMBLY AND START UP PROCEDURE

The bed ships dissasembled mounted on a storage aid.

The assembly is taken care of at the destination by qualified personnel supplied by the operator.

Remove all packaging material and cable ties.

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
ASSEMBLY

THE END BOARDS

The two end boards [3]+[10] must not be interchanged! The adjustment to a lower foot position can otherwise not be achieved.


Look for the decals marking the head and foot sections. They are located in the middle on the crosspipe, in the area of the drive motor mount and in the middle of the the crosspipe on the reclining surface frame [12].

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
ASSEMBLY OF THE RECLINING SURFACE


- Framove the side rail beams [7] and the trapeze pole [2] from the storage aid and put them aside for now.
- Remove the two halves of the reclining surface frame [12] from the storage aid.
- Tut the head side half of the reclining surface frame vertically on the floor. The two trapeze pole sockets [16] point downward while the 2 pairs of knurled screws are pointing upward.
- Now take the foot side half of the reclining surface frame and lift it over the head side half of the reclining surface frame. Now insert the halves of the reclining surface frame into each other.

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- Tighten the 4 knurled screws. Do not use pliers to tighten the screws.
- © Connect the lift rod of the drive motors with the entry of the head side of the reclining surface frame. For that, insert the cotter pins through the receiver and the lift rod and secure it by turning over the clamping yoke.
- Place the combined reclining surface frame flat on the floor.

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
ASSEMBLY OF THE END BOARDS

- Join the end board, head section [10] with the reclining surface frame [12] together. Ensure that the decals correspond!

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Lift the reclining surface frame at the head end and push the two support brackets from the end board, head section [10] into the tubes of the reclining surface frame [12] to the limit.

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There should only be a 5 mm gap between the reclining surface frame and the corner post of the end board.

- Tighten the 4 knurled screws tight. Do not use pliers to tighten the screws.
- Repeat the process with the end board, foot section [3].

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
At the end boards [3]+[10] you will find on the left side and right side a guide bar [17] for the side rail beams [7]. In them runs a plastic slide with 4 metal pins. The metal pins hook into the oblong holes on the face of the 4 side rail

oblong holes on the face of the 4 side rail beams. The side rail beams can be raised or lowered, as needed.

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Draw the two plastic parts of the plastic slide through the long metal strip. To avoid clamping of the plastic slide in the guide track, the metal strip must be positioned precisely around all the

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
parts.

- Position the metal strip [a] above the locking lever over the plastic nose [b].
- Position the lower end of the metal strip [a] directly around the third metal pin [c].
- To mount the side rail beams, the reclining surface frame must be set to the middle position (Electrical connections: see page 17).

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Begin at the end board, head section.

Loosen the knurled screw [d] on the bottom end of the guide track [17] and insert the plastic slide from below.

The release lever [4] must be pointing up.

In order to prevent the plastic slide from falling out of the guide track, reinsert and tighten the knurled screw.

Insert another plastic slide into the opposite side of the same end board.

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
- Repeat the assembly as with the first plastic slide
- Finally, pull the top end of each plastic slide up, at the *head section*, until it locks into place.
- Now place the 8 cm wide side rail beam [7] onto the top pair of metal pins on the end board, head section, and let the other end of the side rail beam lie on the floor.

The side rail beams are rounded on the top. The rounded side must always be uppermost.

- Now place the 7 cm wide side rail beam onto the bottom pair of metal pins and let the other end of the side rail beams lie on the floor.
- Repeat the same procedure on the other side of the end board, head section.

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
Now move to the end board, foot section:

To the end board, foot section, loosen the knurled screws [d] on the bottom end of the guide track [17] and insert both plastic slides from below.

The release lever [4] must be pointing up.

As the plastic slide is being inserted, put in the same time the top pair of the metal pins into the upper side rail beam.

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- At the same time continue to insert the plastic slide upwards in the end board foot section until it locks into place.
- In order to prevent the plastic slide from falling out, reinsert and tighten the knurled screw [d].
- $\ensuremath{\mathscr{F}}$ Repeat the process on the other side of the end board, foot section.
- Check the distance between the beam and the guide track. It should not be greater than 10 mm.

In case the distance is greater, make sure that the support bracket on the end boards assembly is inserted fully into the tube of the reclining surface frame.

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The 4 drive motors are supplied with electricity by the external transformer unit.

Two plugs at the end of the spiral cables must be plugged into the correct drive motors of the end boards.

Before connecting the cables, remove all packaging material.

Take care to correctly connect the respective drive motors for the height adjustment. These must not be interchanged. There should only be one foot level adjustment possible.

Do not place the transformer unit underneath the bed.

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Secure the connecting socket at the end board, head section, in the following manner:

- Guide the connecting socket with the small side into the flap of the strain relief and turn it clockwise until it is horizontally secure.
- Plug in the extension for the drive motor at the **head section** and clip in the strain relief.
- Lay the cables for the drive motors and the hand switch underneath the reclining surface frame
- Plug in the right angle plug for the drive motor of the **foot section** and clip into the strain relief

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Take care that no cable is damaged, that there are no loops, and that the cables are not squeezed by moving parts.

Special care is to be taken when laying the power cable. Make sure it is not damages and there are no loops! The power cable must not be run over by the castors when moving the bed.

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LOCATION REQUIREMENTS

- There must be enough space available to be able to adjust the bed. No furniture, window sills and the like should be in the way.
- The space under the bed must remain clear.
- Check before use on parquet flooring, that the castors do not cause discoloration of the varnish. Use on ceramic tiles, carpet, laminate and linoleum is harmless.

BURMEIER is not liable for any damage to floors due to daily operation.

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- A properly installed 230 volt power outlet must be available, near the bed.
- Take care when attaching further accessory equipment (e.g. compressors for positioning systems, etc.) that their attachment is safe and their functionality is not impaired. Especially important is hereby also the safe placement of all movable connecting cables, tubes, etc. If you have any further questions or need clarification contact the manufacturer of the accessory equipment or BURMEIER.

Observe the following points when setting up the bed to minimize possible fire hazard. Instruct the user and the authorized personnel in these points!

- If possible use only fire-retardant mattresses and bed sheets.
- Only use appropriate firm mattresses according to DIN 13014. This mattress should be inflammable according to DIN 597 Section 1
- Avoid smoking in bed because of the danger of fire.
- Use only technically faultless appliances (e.g. heating blankets) and other electronic devices (e.g. lamps, radios).
- Take care that these devices are only used for the purpose intended and that they cannot accidentally get on top of or underneath the bed sheets (Fire hazard because of overheating).
- Absolutely avoid the use of extension cords or multi outlet extension cords under the bed (Leaking liquids can cause a fire hazard).
- Connect the bed directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.

- REACHING OPERATIONAL READINESS
- After successful assembly, complete an inspection according to chapter 6.
- Clean and disinfect the bed before first time use and before every further use according to chapter 5.

After successfully completing and heeding all steps on page 13 through 19 the bed is now operationally ready.

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
DISASSEMBLY

- Remove the trapeze pole.
- The Lower the bed into the middle horizontal position.
- Lock all four castors.
- Provided in the second of the
- Insert the plastic slide back into the guide bar. Screw the knurled screws tightly back into the guide bars.
- Lower the bed into the lowest horizontal position.
- Pull the power plug out of the socket.
- Unplug the transformer cable from the connecting socket.
- Remove the plugs and strain relief from the drive motors at the head and foot ends.

To disassemble the reclining surface frame proceed as follows:

- Lean the reclining surface frame upright on a wall. Head section downward.
- Remove the cotter pin from the lifting rod of the drive motor for the backrest.
- Loosen the knurled screws on both joints of the reclining surface frame.
- Pull the frame apart.
- Screw the loosened knurled screws back into the reclining surface frame so they do not get lost.

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STORAGE AID


The storage aid joins both headboards and takes in both halves of the reclining surface frame. Additionally there are holders for the side rail beams and the trapeze.

Screw both pieces of the storage aid onto the receiver pin of one of the end boards.


At the same time all knurled screws in the storage aid must face in one direction.

The receptacle for the reclining surface frame must point upward; the basket for the side rail beams go inward and the receptacles for the trapeze point toward the outside.

Use the knurled screws from the reclining surface.

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


- Screw on the second headboard.
- Put the first half of the reclining surface frame (foot section) on the longer receptacle so that the foot section points downward. The mattress brackets [14] point outward.
- Now put the second half of the reclining surface frame (back rest) onto the receptacle, so that the head section points downward. The mattress brackets point inward.

Safely secure the lifting rods of the drive motors with cable ties or something similar.

- Tighten the knurled screws.
- Put the side rail beams into the basket.
- Place the trapeze pole in the appropriate socket.
- The bed is now ready for transport or storage.

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OPERATIONS

GENERAL SAFETY INSTRUCTIONS FOR OPERATION

- In order to avoid toe injuries, wear closed shoes (and not sandals) when operating the bed.
- Make sure that at least three of the castors have their brakes applied.

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SPECIAL SAFETY INSTRUCTIONS FOR THE DRIVE UNIT

- Take care when adjusting the bed that no limbs of patients, users or any other persons, especially children, can be jammed or injured by getting under the movable parts or the bed frame.
- For the safety of the patient, lock the hand switch operation of adjustments to the back and thigh rests to avoid unwanted adjustments when the side rail beams are up. (Danger of getting limbs trapped when adjusting the back- and thigh rests).
- When using accessories with the bed it is necessary to: make sure that the arrangement of these accessories will not cause any places that might squeeze or hurt the patient. If this is not the case the operator must safely disable the adjustment function of the back- and thigh rests.
- ▶ To do so place the hand switch out of the patient's reach (e.g. foot section) or lock the adjustment functions on the hand switch.
- Make sure that the power and hand switch cable cannot be jammed or damaged in other ways.
- Ensure that only hand switches that are in perfect condition are used. Faulty hand switches (e.g. with torn or indented foil) can result in spontaneous operation of the bed and can thus endanger the patient.

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- Make sure that before relocating the bed, the power cable cannot be stretched, rolled over or damaged in any other way.
- Before every relocation, unplug the transformer unit cable from the wall.
- Position the transformer unit in a safe place so that it cannot fall off.
- Ensure that the connecting socket is always closed with the stopper when it is disconnected from the transformer unit.
- Make sure that when adjusting the bed no obstacles such as furniture or slanted ceilings are in the way. This will prevent damage.

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- ◆ An electronic overload detector switches off the central control unit, when the load is too great. After clearing the overload and pressing the hand switch again the central control unit will work again.
- A continuous run time of 2 minutes should not be exceeded. Afterwards a pause of no less than 18 minutes should be observed. Another possibility: 1 minute continuous run time and 9 minutes pause.
- If the maximum continuous run time is not heeded (e.g. by constantly "playing" with the hand switch) a thermal safety component switches off the power for safety reasons.
- The adjustable range of every function is electrically/mechanically possible only within a predetermined range.
- Electrical disturbances affecting other nearby electronic appliances (e.g. crackling on the radio) cannot be excluded, despite adherence to all prescribed regulations. In such rare circumstances increase the distance between the appliances, do not use the same power socket or switch the affected appliance off for a while.

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HAND SWITCH

The electrical adjustments can be made using the hand switch [8].


The adjustable range of every function is electrically/mechanically possible only within a predetermined range.

Locking mechanisms have been built into the hand switch for safety reasons. For the safety of the patient, adjustment functions can be locked at the hand switch, if, from a clinical point of view and with the approval of the attending physician, the condition of the patient warrants

- ◆ The hand switch, with its elastic hook can be hooked onto the bed.
- The spiral cable allows for greater mobility.
- ◆ The hand switch is water-tight and washable (IP X4).

OPERATING THE HAND SWITCH

- Train the patient in the functions of the hand
- The drive motors will run as long as the button is held down.
- The adjustment functions, except for the lowered foot position, can be used in both directions.
- Basically the functions of the buttons are as shown:
 - raise
 - lower

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Adjusting the backrest

Push this button if you want to adjust the backrest angle.

Adjusting the reclining surface height

Push this button if you want to adjust the height of the reclining surface.

If the bed is in the lowered foot position, it will automatically be adjusted to horizontal level, when it is readjusted to the lowest or highest position.

Adjusting the thigh rest

Push this button to adjust the thigh rest angle.

Adjusting to the lowered foot position

Push this button if you want to adjust to the lowered foot position. With this button only the function "Down" is possible.

If the bed is in the lowered foot position, it will automatically be adjusted to horizontal level, when it is readjusted to the lowest or highest position.

■ Unlock the foot brake of both castors on the head or foot section to avoid damaging the floor.

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LOCKING FUNCTION

Only the operator is authorized to use the locking function!

In case the clinical condition of the patient is so critical that any adjustment of the beds by means of hand switch causes danger, the operator has to disable the bed immediately. The nursing bed remains in the position it had at the time of disconnecting.

Do not turn the locking key with force over the stop of the lock! The lock or the whole hand switch can be damaged.

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The locking key is delivered attached to the bed with a cable tie on the hand switch. The locking key is not intended for use by the patient. It must be removed from the hand switch. The operator or a person authorized by the physician should keep the **locking key**.

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■ Turn the respective lock on the hand switch with the locking key clockwise in the locking position.

The color of the respective display will change from green to yellow.

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Lock is vertical

Color of the display is: green

Buttons can be pushed ("Clicking sound")

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
Drive disabled

Lock has been turned clockwise by about 15°

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CASTORS

The bed stands on four steerable, individually-locking castors [11].

- Braking: Push the foot lever down.
- **Driving**: Push the foot lever toward the corner pole.

- Tragen Sie beim Bedienen des Bettes geschlossene Schuhe, um Verletzungen an den Zehen zu vermeiden.
- Stellen Sie sicher, dass mindestens drei Laufrollen gebremst sind.
- In order to avoid toe injuries, wear closed shoes (and not sandals) when operating the bed.
- Make sure that at least three of the castors have their brakes applied.

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TRAPEZE POLE

■ The maximum load of the trapeze pole is 75 kg at the front end.

The trapeze pole attached to the bed makes it easier for the patient to get into and out of bed.

In both corners of the head section of the reclining surface frame are two trapeze pole sockets [16]. The trapeze pole socket has a notch at the top [c], which together with the metal pin [b] restricts the swivel range of the trapeze pole. The trapeze pole should be installed on the side where the patient gets in and out of bed.

Attaching

Insert the trapeze pole [2] into the socket.
The metal pin [b] must be in the notch [c] of the socket.


Removal

Pull the trapeze pole straight up out of the socket.

Trapeze pole swivel range

Do not swivel the trapeze pole outside of the bed. There is a danger that through the pull, the bed will tip over.

The metal pin must always be in the notch of the socket.

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TRAPEZE

A trapeze (grab handle) can be attached to the trapeze pole. The patient can hold on to the trapeze to make it easier to change to a different position.

Regularly check the trapeze at the strap for any damage (see Chapter 6). A damaged trapeze or strap must immediately be exchanged.


Durability

There is a date clock on the trapeze. With normal use the trapeze will last for at least five years. Afterwards a sight check is to be performed periodically to determine if further use is permissible.

Adjustable range of the trapeze

The height of the trapeze can be adjusted by means of the strap to accomodate a range of about 55 cm to 70 cm (measuring from the top of the mattress).

- Push the loop of the trapeze [1] over the first stud of the patient grab handle pole.
- Test by tightly pulling the trapeze downward that it is secure.

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Useful Tip: The maximum lifting load of the trapeze pole is 75 kg at the front end.

- The height of the trapeze can be adjusted by means of the strap.
- Make sure the the strap is fed correctly through the buckle.
- Make sure that the end of the strap is at least 3 cm above the buckle.

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LOWER LEG REST (OPTIONAL)

When the thigh rest is fixed, the lower leg rest [5] can be individually adjusted. For this reason there are two adjusting ratchets, which will lock the lower leg rest into position.

Raising

Raise the lower leg rest on the foot end - not with the mattress brackets - up to the desired level. The lower leg rest will automatically snap into place.

Lowering

- First lift the lower leg rest all the way to the stopping point.
- Slowly lower the lower leg rest.

There is a potential for injury if the lower leg rest falls straight down.

If you lower the thigh rest, the lower leg rest will be lowered automatically.

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SIDE RAILS


Side rails conveniently prevent the patient from falling out of bed accidentally. They do not, however, prevent a patient from leaving the bed.

Engaging

- Pull the side rail beams [3] on one end after the other until until they lock in its highest position. You should not be able to move them up or down.
- Check that the side rail beams are locked into place by pushing on them from above.

Disengaging


- Lift the side rail beams slightly.
- Push the upper release lever [4] and slowly lower the side rail beams.

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Special safety notes for use of side rail

Only use safely lockable, technically perfect, undamaged side rail beams with proper gap dimensions.

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- Make certain that only suitable Original BURMEIER side rail beams are used, which have been approved for use with the bed model you are using or that are already integrated in the bed.
- Before using side rail beams evaluate and take into account the clinical condition and any special problems of the patient
- ▶ If, for example, the patient is very bewildered or very anxious, forgo the side rails and use alternative safety measures such as a restraining belt or restraint sheets.
- ▶ For especially small, frail patients, additional protective measures for reducing the distance between bars on the side rails may be necessary. In these cases use, for example side rail foam covers (accessories, see chapter 8), retaining belts, etc.

This will effectively ensure patient safety by reducing the risk of becoming trapped or slipping through.

- Only use appropriate, firm mattresses accord to DIN 13014 with a height of 10 to a maximum of 18 cm. Further, this mattress should be fire resistant according to DIN 597 Section 1 and 2.
- If elevated patient positioning systems are employed, such as antidecubitus mattresses, an effective side rail height of at least 22 cm above the unweighted mattresses must be ensured. If this requirement is not met, the operator must take suitable measures to safeguard the patient. This can be done by using an attachable bar for side rails (Accessory).
- The mechanical parts on the bed frame and the side rails which serve to protect the patient must be regularly examined for damage. This should not only happen before the installation of the side rail but also during use in shorter intervals (at least before every new assignment).
- To safeguard the patient against unintentional adjustments, always lock the electrical adjustment of the back and thigh rests on electrical beds when the side rails are raised. Otherwise there is a risk of trapping limbs when the back and thigh rest are adjusted.
- Before lowering the side rail beams make sure that the patient's limbs are clear of the side rail beams. This way you can avoid injuring the patient.

If these safety warnings are not observed there is a danger for the patient of being injured by the side rails, of falling out of bed because of too large a distance between the bars or because of damage, the use of unsuitable side rails or the incomplete locking of the side rail.

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EMERGENCY LOWERING OF THE BACK REST

The upright back rest can be lowered manually in case of a power outage or failure of the power supply or the electrical drive systems.

In this case two users are absolutely necessary!

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

Disregarding these safety guidelines and instructions may cause the backrest to fall out of control, which could lead to serious injuries for the user and the patient!

Emergency lowering may only be carried out in the case of extreme emergency and by users who have a complete command of the procedures below.

We strongly advise you to practice emergency lowering of the backrest under normal conditions so that in case of an emergency you will be able to react quickly and properly.

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If the back rest falls out of control, the patient and the second user may be injured!

To convert the bed back to its normal operating state: Swing up the connecting mount and secure the spring pin, then turn over the bent bracket.

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CLEANING AND DISINFECTING

GENERAL

Cleaning is the most important requirement for a successful chemical disinfection.

Generally a routine cleaning of the bed, if used by the same patient, is hygienically sufficient.

Disinfection of the bed frame is only necessary in case of visible contamination with infectious or potentially infectious material (blood, stool, pus) or by presence of an infectious disease under the direction of a physician.

The bed must be cleaned and washed down with disinfectant in case of a patient change!

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Before beginning the cleaning, please note:

- Unplug the power cable from the wall socket and store it in such a way that it does not come into excessive contact with water or other cleaning solutions (store in plastic bag).
- Ensure that all plugs on the transformer unit and the drive motors are plugged in properly.
- None of the electrical components must show any visible damage; otherwise water or cleaning solution could penetrate. This could lead to malfunction or damage to the electrical components.
- Before operating the bed again, make sure that there is no residual moisture on the electrical contacts, by drying or blowing on the power plug.
- The electrical components must not be cleaned with a water jet, a high pressure cleaner or other similar device! Clean only with a moist towel.
- If you suspect that water or any other form of moisture has penetrated the electrical components, pull the power plug out of the wall socket immediately. If already disconnected from the power supply, make sure it is not plugged in again. Clearly mark the nursing bed as "Defective" and take it immediately out of service. Report this to the operator promptly.
- Failure to follow these directions could lead to major damage of the bed and its electrical components and lead to subsequent malfunctions.

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CLEANING AND DISINFECTION PI AN

- Remove the bed sheets and give them to the laundry.
- Clean all surfaces including the slatted boards and the reclining surface made of plastic inserts or wire mesh base with a mild and environmentally safe cleaning solution. This also applies to the hand switch.
- If the bed has visible contamination e.g. infectious or potentially infectious material has polluted it, it must be disinfected by wiping after cleaning. Use approved disinfectants for each surface according to the products listed by the DGHM (German Organization for Hygiene and Microbiology). This also applies to all beds from patients with illnesses subject to report according to §6 of the infection control law (German abbreviation: IfSG, Protection against Infection Act), colonization or infection with treatment resistant pathogens (e.g. MRSA, VRE) and all beds from the intensive care and infection control units. Use only the DGHM approved disinfectant in the approved concentration.
- Disinfecting the castors is only necessary when visibly contaminated by infectious or potentially infectious material.

Notice: Continuous disinfection is only necessary for patients with a treatment resistant pathogen (e.g. MRSA) within the hospital.

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INSTRUCTING THE USER AND QUALIFIED PERSONNEL

To guarantee the correct procedure for cleaning and disinfection is followed, we recommend to correctly instruct the handler and qualified personnel. Thereby communicating that the following points need to be observed:

- The clean bed is to be transported to the patient's home in such a way that it does not become polluted or contaminated.
- When disassembling the bed we recommend cleaning and disinfecting it right away. The qualified personnel should be familiar with the steps necessary to thoroughly clean and disinfect (the operator will provide the individual steps necessary). It is necessary to use only the DGHM approved disinfectant in the approved concentration. The disinfectant must be appropriate for the surface being cleaned.
- The qualified personnel must be equipped with water-tight one-time-use aprons and gloves.
- Only fresh and clean towels must be used to clean and disinfect, these must then be given to the laundry.

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■ After cleaning the bed, the qualified personnel must also disinfect their hands before going to the next patient. An appropriate disinfectant dispenser should be part of the qualified personnel's gear.

The immediate cleaning of the bed while disassembling has the benefit of not transporting "dirty" beds or bed parts along with clean ones in the same truck. A transfer of potentially infectious spores which may be clinging to the used bed may be avoided. A transfer of spores of a hospital-acquired infection will be avoided with strict adherence to these recommendations.

If the beds are to be stored for a while, they should be covered to keep them free of dust, accidental pollution and contamination.

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CLEANING SUPPLIES AND DISINFECTANTS

Make use of these recommendations to keep the bed serviceable for a long time.

- The surfaces must be intact; any damage should be repaired immediately.
- We recommend a (damp) wash down. When choosing the cleaning solution make sure it is mild (gentle to skin and surfaces) und environmentally friendly. In general a regular household cleaner can be used.
- For cleaning the coated metal parts use a wet cloth with regular household cleaners.
- The disinfectant used for the wash down should be one approved for the purpose and listed by the DGHM (in the recommended dosage). Disinfectants and concentration listed by the Robert-Koch-Institute (the former Bundesgesundheitsamtes) should only be used in case of an epidemic according to the infection safety law and only on the orders of the appropriate health department.
- Despite the excellent mechanical resistance of the bed, scratches, markings, etc. which penetrate the coating should be resealed using a suitable medium to prevent the penetration of moisture. Contact BURMEIER or your specialized dealer with any questions.

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■ Do not use scouring agents, stainless steel cleaners, abrasive cleaning agents or scouring pads. These substances can damage the surface.

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Notice:

Surface disinfectants based on aldehyde generally have the advantage of working on a large array of pathogens, have a relatively low protein error and are environmentally friendly. Its main disadvantage is the allergy-causing and irritation potential. Compounds on the basis of glucoprotamine do not have this disadvantage and are very effective, though mostly more expensive. Disinfectants on the basis of chlorine splitting combinations, with prolonged contact or in too high doses, can be corrosive to metals, plastics, rubber and other materials. In addition to this, the compounds have a high protein error, are irritants to the mucous membranes and are not environmentally friendly.

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CONTACT WITH DISINFECTANTS

- Follow the guidelines for each product and manufacturer. Observe the exact dosage! We recommend automatic dosage apparatus.
- Always mix the solution with cold water, to avoid vapors that aggravate the mucous membranes.
- Wear gloves to avoid direct skin contact.
- Do not store the mixed disinfectant in open containers with the cleaning rag inside. Close the container!
- Use closeable bottles with dosage pumps for wetting the cleaning rag.
- After finishing the surface disinfection ventilate the room.
- Use a wipe down disinfecting method, do not use a spray. Using the disinfectant as a spray causes the disinfectant to turn into a mist which could be inhaled.
- The wipe-down disinfection plays a major part.
- Do not use alcohol for larger areas.

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operators of nursing beds are obligated to keep medical products in a safe state throughout their entire service life. This also includes regularly carrying out expert maintenance as well as safety checks.

If a bed is purchased for domestic use (not for commercial use), the regular safety inspection is not compulsory but is still recommended by the manufacturer.

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Notice for Operators

This bed has been constructed and built to last safely for a long time. With appropriate operation and use, this bed has a life expectancy of 2 to 8 years. The life expectancy is based on the conditions and frequency of use.

Through repeated transport, assembly and disassembly, improper operation as well as long time use, it is possible that damage, defects and deterioration can occur. These deficiencies can lead to danger if they are not recognized and quickly fixed.

Because of this there are legal bases for the regular enforcement of tests to ensure the safe condition of this medical product.

According to § 4 of the Medizinprodukte-Betreiberverordnung operators are obliged to maintain the product. Because of this the following regular inspection and function controls are to be performed by the operator as well as the user.

Point out these necessary tests to the user.

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BY THE USER

Along with the regular substantial tests by technically qualified personnel, the normal user (nurses, supporting relatives, etc.) must also perform a sight and function check at shorter, regular intervals, as well as before a new patient.

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Recommendation:

All electrical and mechanical components should be tested monthly. The main power cable as well as the hand switch cable should be checked after every mechanical load and after every relocation.

For this purpose use the checklist on page 35.

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Checklist: Test protocol by the user

Test		ОК	Not OK	Defect description	
Visual check of electrical components					
Hand switch and cable	Damage, Cable Routing				
Cable harness	Damage, Cable Routing				
Hand switch	Damage, Foil				
Transformer unit	Damage, Cable Routing				
Visual check of mechanical componen	ts				
Trapeze pole/Socket	Damage, Deformation				
Bed frame	Damage, Deformation				
Spring wood slats	Damage, Splinters				
Bed box	Damage, Splinters				
Reclining surface frame	Damage, Deformation				
Side rail	Damage, Splinters				
Function test of electrical components	Function test of electrical components				
Hand switch	Functiontest, Blocking function				
Function test of mechanical components					
Emergency lowering of backrest	Test according to instruction manual				
Castors	Free movement an locking				
Knurled screws	Tight fit				
Side rail	Safe locking/unlocking				
Lower leg rest	Latching				
Accessories (e.g. trapeze pole, trapeze)	Secure, Damage				
Signature of the tester:	Result of the test:			Date:	

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BY THE OPERATOR

The operator of the nursing bed according to MPBetreibV § 4 is obligated with every new assembly, every maintenance and throughout the operation, to perform regular inspections, to ensure the safe condition of the nursing bed.


These tests are to be repeated during regular maintenance and depending on the use according to MPBetreibV § 4 and from the Directive of the German Employers Liability Insurance Association according to BGV A3.

- ◆ Observe this sequence for the inspection according to DIN EN 62353:
- I. Visual inspection
- II. Electric Measurements
- III. Function Test
- ◆ The performance of the function test as well as the rating and documentation of the test results can, according to MPBetreibV § 4 only be performed by experts with the appropriate knowledge, training and experience as well as the tools necessary to perform it according to the rules.
- ◆ The performance of electrical measurements according to DIN EN 62353 can, with availability of proper measuring equipment, also be performed by an electro-technically trained person [according to BGV A3] with medical and appliance specific training.
- ◆ The assessment and documentation of the test results can only be performed by an expert with the corresponding knowledge, training and experience.
- The testing of the transformer unit only requires measurement of the leakage current.
 The bed can therefore be used again once the transformer has been replaced with a new unit.
- BURMEIER offers transformer testing as a customer service. Send the transformer unit back to BURMEIER and receive a tested replacement unit. Contact us for further details; Address see page 42.

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Preperation:

- Unplug the power plug of the transformer unit from the wall socket.
- Unplug the 24 volt cable of the transformer unit from the connecting socket.
- Plug the plug of the 24 volt cable in the adapter (extra accessory; request from BURMEIER).
- Connect the adapter with the test probe of the test unit.
- Plug the transformer power plug into the power socket of the test unit.

Test method:

- Output power test: direct or differential according to DIN EN 62353
- Output power test according to the test unit manufacturers instructions.

Limit:

■ Output power less than 0.1mA.

Test cycle:

We recommend a yearly visual inspection and performance check. If this checks was been ok, an electrical testing suffices every ten years.

If the 2% error quota (see BGV A3: §5, table 1 B) is shown to have been observed, the test cycle for inspections and functional testing can be extended to a maximum of two years.

If damage or a function error is suspected take the slatted frame out of service immediately and disconnect it from the power until an exchange or repair can be made!

The test protocol draft on **pages 40 and 41** can be used.

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Customer / med. Facility / practice:							
Address:							
Carr	ied out: ☐ Repeat inspection			☐ Inspection	n prior initial operation		
				☐ Inspection	on following repairs / servicing		
Equipment type: ☐ Hospital bed ☑ Nursing bed				Protection classification: □ ☑			
Bed	type: DALI LOW ENTRY 24 Volt			Inventory Number:			
Loca	ntion:			Serial Number:			
Grad	le level (IEC60601-2-52): 1 1 2 3 4 5						
Manufacturer: Burmeier GmbH & Co. KG			User-specific: none				
Testing equipment used (Type / Inventory Number):			1.				
MPG Classification: Class /			2.				
I. Visual Inspection				ок	Not OK	Defect Description	
Visu	al inspection of the electriacl components						
Stick	ers and type plates	Presence, legib	le				
Tran	sformer unit housing	Damage					
Moto	r housing and lifting rods of motors	Correct position	n, Damage				
Moto	r bolt	Check that the	lock washer is pr	operly in place			
Han	d switch	Damage, Foil in	ntact				
Cabl	e for motor and hand control	Damage, Routing					
Plug	of harness / switching power supply	Available, Correct position					
Visu	al inspection of the mechanical components						
Stickers and type plates Presence, legible							
Safe	fe working load / patient weight Presence, legible						
Trap	Trapeze pole / socket Damage, Deformation						
Bed frame Damage, Deformation							
Spring wood slats Damage, Splinters							
Castors Damage							
Reclining surface Damage, Deformation							
Woo	Wooden boards Damage, Splinters						
Weld	Welded seams Split welded seams						
Side	e rail Damage, Splinters						
Knur	inurled screws Tight fit						
Wea	Wearing parts such as joints Damage						
Lock	ing ring of backrest/thigh rest	Check that the	ring is properly in	place			
II. Electric measurments according to DIN EN 62353:2008-08 (leakage current, direct)							
Proceed as follow:				ОК	Not OK	Defect Description	
1.	Plug in transformer unit in test socket of testing device						
Insert the plug of the 24 volt cable into the special adapter (extra accessory; request from BURMEIER)							
Connect the special adapter with the test probe of the test unit.							
4. Start procedures with the measuring instrument							
	Limit Value Value of first measurement value			Current actual value			
Result: Bed SK II (Type B) 0,1 mA mA			mA				


Test protocol of an inspection of electro-medical equipment according to DIN EN 62353: 2008-08 Page 1 of 2

Not III. Performance check OK **Defect Description** ОК Performance check of the electrial components Motion limit switching of the motors Automatic cut-off Hand switch, blocking functions Test according to instructions for use Drive motors Abnormal Noise Development Drive motors Test according to instructions for use Strain relief for power cable Power cable firmly fastened Performance check of the mechanical components Joints and pivots Smooth operation Trapeze with belt strap Secure hold Castors Drive and brake Emergency lowering of the backrest Test according to instructions for use Side rail Safe locking / Secure hold / Safe unlocking Lower leg rest Latching Accessories (e. g. trapeze pole, trapeze) Secure, Damage, Suitability Result of the inspection: All values within permissible range: ☐ yes □ no Next inspection date: Inspection passed: □ yes \square Defective, do not use the bed! => Repair If inspection was not passed: ☐ Defective, do not use the bed! => Take out of Service ☐ Bed does not meet the safety standards Test approval sticker applied: ☐ yes □ no Part of this inspection report: Comments: Inspected on: Inspector: Signature: Evaluated on: Operator / Expert: Signature:

Test protocol of an inspection of electro-medical equipment according to DIN EN 62353: 2008-08 Page 2 of 2

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SPARE PARTS

The corresponding replacement parts can be obtained from BURMEIER by specifying the item number, the order number and the serial number. The necessary details are found on the type plate under the cross tubing of the reclining surface frame.

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MANUFACTURER'S ADDRESS

In order to maintain operational safety and the right to claim under warranty, only original BURMEIER replacement parts may be used!

For ordering replacement parts, customer service requests and other questions, please contact:

Burmeier GmbH & Co. KG (A subsidiary of the Stiegelmeyer Group)

Pivitsheider Strasse 270 32791 Lage/Lippe, Germany Phone: ++49 52 32 / 98 41-0 Telefax: ++49 52 32 / 98 41-41

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REPLACEMENT OF ELECTRIAL COMPONENTS

Dangerous voltage, electric shock possible

- Before doing any work unplug the main power cable from the wall socket!
- Any work on the electrical equipment may only be performed by customer service, the motor manufacturer or qualified and authorized electricians under consideration of the applicable VDE rules and safety guidelines!

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Exchange of the Hand switch

- Unplug the power plug from the wall socket.
- Unplug the plug at the end of the spiral cable of the hand switch from the connector underneath the reclining surface.
- Lay the cable of the new hand switch so that it cannot be damaged by moving parts of the bed.
- Insert the plug at the end of the spiral cable of the new hand switch into the connector of the hand switch.

Take care that the O-ring on the plug is not damaged, otherwise the connection may not be moisture-proof.

Tarry out a function test of the electrical adjustments after the exchange.

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TROUBLESHOOTING

The following table on page 45 is a guide to rectifying common malfunctions:

If a function fault is detected during operation, which the user cannot rectify with the use of the trouble-shooting table aid, the specialized staff for maintenance and repair of the appropriate owner is to be called.

The user may **under no circumstances** try to eliminate defects at the electrical equipment!

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Dangerous voltage, electric shock possible!

Any work on the electrical equipment may only be performed by service personnel, the motor manufacturer or qualified and authorized electricians as required by the applicable VDE rules and safety guidelines!

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Trouble-shooting table aid

Problem	Possible causes	Solution			
Hand switch / drive system not functioning	■ Hand switch plug not properly connected				
	■ Transformer plug (24 volt) not inserted properly	Press the transformer unit plug firmly into the connecting socket			
	■ Hand switch or drive system defective	Inform your operator about necessary repairs			
	■ Functions are blocked on hand switch	Unlock functions (see page 23)			
Green LED on transformer does not light up Power cable not properly connected		Plug in the power cable properly			
and drive systems do not function	■ Wall socket without power	Test the wall socket and fuse box			
	■ Transformer defective	Exchange the transformer; inform your operator about any necessary repairs			
Green LED on trans-	■ Hand switch is defective	Exchange the hand switch			
former lights up continuously, yellow LED does not light up when any	■ Transformer plug (24 volt) not inserted properly	Test the plug connection			
hand switch button is pressed	■ Functions are locked on hand switch	Unlock the functions (see page 23)			
Drives operate in one direction only	■ Hand switch or drive system is defective	Inform your operator about necessary repairs			
Drives stop suddenly after prolonged adjust-	■ Heat swich in transformer is activated	Pay attention minimum pause (see page 22)			
ment times		Exchange the transformer; inform your operator about any necessary repairs			

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ACCESSORIES

Only the Original BURMEIER accessories may be used. If other accessories are used, we can not be liable for any resulting injuries, defects and dangers.

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Make sure that the arrangement of accessories does not produce any crush or shearing zones for the patient when the back and leg rests are adjusted. If this cannot be ensured, the user must safely prevent the patient from adjusting the back- and leg rest.

▶ To do so, attach the hand switch out of reach (e.g. at the foot end of the bed), or lock hand switch adjustments.

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Disassembled bed:

Chassis complete incl. motor

Reclining surface frame with motors

4 side rail bars, wooden 4 side rail bars, iron

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SPECIAL FEATURES

- Electrical height adjustment of the reclining surface from about 40 cm to 80 cm.
- Electrical adjustment of the backrest from 0° to about 70°.
- Electrical adjustment of the leg rest from 0° to about 35°.
- Electrical inclination of the reclining surface to the foot end from about 10°.
- Driven with four individually-locking castors.
- ◆ Reclining surface of 200 x 90 cm, divided into four areas; Outside dimensions are about 212 x 101 cm.
- Side rails on both sides can be lowered.

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Transformer unit

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Input Voltage Max Amps

Output Volts
Output Amps
Operating time

Classification, Protection class II, Type B, MPG classification Class I, not for use in explosive atmospheres.

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Linak HL 74

IP X4 Schutzart

Typ

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Linak LA 24

1500 N / 575 mm / 405 mm

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DC 24 V

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> IP X4 Schutzart

> > Motor Rückenlehne

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Linak LA 24

3500 N / 320 mm / 110 mm

Kraft / Einbaumaß / Hub DC 24 V

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> IP X4 Schutzart

> > Motor Oberschenkellehne

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Linak LA 24

2500 N / 270 mm / 60 mm

DC 24 V

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IP X4

Schutzart

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Hand switch with locking function

Type

Protection

Motors for reclining surface height

adjustment

Type

Power / Size / Lift Input Voltage Operating Time

Protection

Motor for backrest

Type

Power / Size / Lift Input Voltage Operating Time Protection

Motor for thigh rest

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Power / Size / Lift Input Voltage Operating Time Protection

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Zulässiges Patientengewicht

Explanation of Icons:

 Unit with Type B application according to IEC 601-1

(Special protection against electric shock)

- ◆ Protection Class II device, double insulation
- Transformer with thermal fuse block
- Safety transformer in accordance with VDE 05510
- ◆ Caution! Consult instruction manual
- ◆ The electrical equipment offers protection from internal dust deposits and is splash-proof from all sides
- The electrical equipment is dust-proof from internal dust deposits and is splash-proof from
- Conformity mark in accordance with Medical Products directive 93/42 EWG
- · Safety weight load
- Permitted patient weight

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Relative humidity

Air pressure

Noise level during motor adjustments

The ambient conditions stated below must be

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IEC 60601-1-2-52: 2009-12 Medical electriacl equipment:

> Particular requirements for the basic safety and essential performance of medical beds

GUIDELINES

Risk analysis for medical devices

Safety for electromedical equipment

Electromagnetic compatibility Medical electrical equipment:

with fitness for use

Medical electrical equipment – Requirements for medical electrical equipment and medical electrical systems used in the home healthcare

environment

Particular requirements for the safety and essential performance of medical beds

Medical electriacl equipment:

Particular requirements for the basic safety and essential performance of medical beds

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DISPOSAL INSTRUCTIONS

- This nursing bed is if electrically adjustable classified as a commercially used electric appliance (b2b) in accordance with WEEE-Guidelines 2002/96/EG (Electronics Law).
- ◆ Exchanged electronic components (drives, control unit, hand switch, etc.) of these nursing beds are considered electronic scrap and must be disposed of according to the WEEE-Guidelines.
- With nursing beds that came into use after August 13, 2005 (see designation "Serial" on the type label at the head of the bed), the operator is required to not dispose of the electrical components in the municipal waste area, but to send them directly to the manufacturer. BURMEIER and their service and sales partners will take these parts back.
- Our general terms and conditions apply for this redemption.
- No longer useable rechargeable batteries have to be disposed of according to regulations and do not belong in the house garbage.
- With all disposed components the owner must make sure that they are not infectious or contaminated.
- In the event that the nursing bed must be scrapped, the plastic and metal parts used must be disposed of separately, following the proper disposal procedure.
- Pneumatic springs are under high pressure! Before disposal they have to be depressurized. Please contact the manufacturer of the device on how to depressurize them (see type plate).
- Hydraulic units are made up of mostly metal and plastic. Before disposing according to manufacturer guidelines the oil has to be drained and disposed of properly.
- Contact your local community, disposal facility or our service department with any questions.

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Internationale Normen:

• IEC 60601-1-2-52 : 2009-12 Medical electriacl equipment: Particular requirements for the basic safety and essential performance of

medical beds

Lage, den 2015-05-28

(Kampisiulis Kemmler)

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(Rekemeier)

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EC-DECLARATION OF CONFORMITY

We,

Burmeier GmbH & Co. KG Pivitsheider Strasse 270 32791 Lage / Lippe Germany

hereby declare under sole responsibility as the manufacturer that the product model named below:

Nursing bed DALI LOW ENTRY 24 Volt

in the version submitted complies with the regulations of the EC Directive 93/42/EEC for Medical Devices, last amended by Directive 2007/47/EC dated 5 September 2007. It is categorised as a Class I active medical device.

The relefant technical documentation is kept by the manufacturer's safety representative.

To evaluate the conformity to the Directives, all applicable parts of the following standards were referred to:

complies with the regulations of EC Guidelines 93/42/EWG for medical products and fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for care beds.

Harmonised standards:

EN 14971: 2013-04 Risk analysis for medical devices
 EN 60601-1: 2007-07 Safety for electromedical equipment
 EN 60601-1-2: 2007-12 Electromagnetic compatibility

• DIN EN 60601-1-6 : 2010-10 Medical electrical equipment: with fitness for use

DIN EN 60601-1-11: 2010-10 Medical electrical equipment – Requirements for medical electrical equipment and medical electrical systems

used in the home healthcare environment

DIN EN 60601-2-52: 2010-12 Particular requirements for the safety and essential performance of medical beds

International standards:

• IEC 60601-1-2-52 : 2009-12 Medical electriacl equipment: Particular requirements for the basic safety and essential performance of medical beds

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Email

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No part of this manual may be reproduced without the prior written permission of the publisher!

All rights reserved.

Subject to technical changes!

CE

Stand: 2015.08.28 Updated: 2015.08.28

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Burmeier GmbH & Co. KG (Ein Unternehmen der Stiegelmeyer-Gruppe)